

If your laboratory is helping your community return to school, Quiksal Oral Rinse is an ideal solution for student testing.

The Quiksal Oral Rinse Specimen Collection Device is designed to replace nasal swabs and saliva collection devices for use in COVID-19 testing. This collection device is a 2ml oral rinse that is swished for 60 seconds and returned to its container.

- Non-invasive and can be self administered by the student
- Post stabilization buffer extends stability and deactivates the virus
- Safer handling and pre-analytical processes
- Reduced viscosity makes testing easier and reduces cross-contamination
- Food grade ingredients stimulate the salivary glands eliminating dry mouth concerns
- Reduced non-infectious positives
- Optimized for specimen pooling
- Validation specimens available

If you would like more information on how we can validate your laboratory using this method, please reach out to Community Care at 844.419.2470 or purchasing@communityhcps.com

COVID - 19 Oral Rinse Instructions

Step 1



Patient should have nothing by mouth 15 minutes prior to collecting sample.

Step 2




Write the patient's full name, date of birth, and date/time of collection on both labels (top and side) on tube. Please top label with patient info onto the test order form.

Step 3



Remove the tube cap from the collection tube and hand to patient.

Step 4




Instruct patient to dispense the oral rinse fluid from the collection tube into mouth. Return tube to collector. **DO NOT swallow contents.** Instruct the patient to swish the contents around vigorously for a full minute. After 1 minute have patient spit the entire contents into collection tube.

Step 5



Add the stabilization solution from the syringe (**Note: remove cap from syringe first**) or tear-away pouch directly into the collection tube *after* patient has spit in tube.

Step 6



Carefully place the cap securely on the tube. Close tightly to prevent leaking. Invert tube 2-3 times to fully mix stabilization solution with oral rinse from patient.

Step 7



Place collection tube into bio-hazard zip-closure bag and seal completely. Place the folded test order form in the front pouch of the bio-hazard bag.

FAQ:

1. Patient should not take anything in the mouth 15 minutes prior to sample collection. This includes food, gum, any tobacco products, fluids and medications. If the patient has consumed anything by mouth, you must wait 15 minutes after consumption before performing the test.
2. If the patient accidentally swallows the oral rinse, instruct them to spit out the contents and start over. Swallowing the oral rinse does not cause harm, it is completely food grade.
3. The stabilization solutions in step 5 may come in the form of a syringe or a tear-way package. Do **NOT** put stabilization solution in mouth.
4. Testing kits should be stored in cool and dry place at 35 - 86 F degrees.

Thank you for purchasing the COVID-19 Oral Rinse kit from Community Care. We are here to answer any of your questions, please contact us at Community Care Client Services.

CS@elementhcp.com

844-419-3470